



*Laysan Bio Inc.*

## **PRODUCT STORAGE INSTRUCTIONS**

Because PEG reagents are temperature, air and moisture sensitive they should be stored in a freezer at or below -10C upon receipt. Please allow PEG reagents to equal room temperature before opening. We recommend storing our PEG reagents as the dry, powder form only and not in any type of solution. Each time the bottle is opened it will need to be backfilled with an inert gas (i.e. Argon, Nitrogen) and tape the lid before being stored again.

**WE DO NOT RECOMMEND YOU STORE PEG PRODUCTS IN SOLUTION.**

*If your reagent requires rebottling please follow these steps:*

Allow reagent to equilibrate to room temperature

While reagent is thawing prepare glove box or another suitable area for packaging which will allow for an inert atmosphere

Flush the repackaging enclosure with an inert gas (i.e. Argon, Nitrogen)

Repackage into desired aliquots

Flush each aliquot with inert gas before closing and tape the lid

Store each aliquot at or below -10C until ready to use

If you do not have the proper equipment to rebottle Laysan Bio offer a bottling service at \$50 set up fee and \$5.00 a bottle. This service must be ordered before the product ships.

If you have questions regarding this information we encourage you to contact our sales team at 256.586.9004 or [sales@laysanbio.com](mailto:sales@laysanbio.com).

For technical assistance email [tech@laysanbio.com](mailto:tech@laysanbio.com).

**Laysan Bio Inc.**  
**1560 Tower Drive**  
**Arab, AL 35016**  
**Phone: 256.586.9004**  
**[sales@laysanbio.com](mailto:sales@laysanbio.com)**

1560 Tower Drive Arab, Alabama 35016  
Phone: 256.586.9004 Fax: 256.586.9007